

How to Set Up Continuous EEG (C EEG)

OBTAIN SUPPLIES

1. *EEG module (yellow)*
2. EEG cable with wires
3. NuPrep™ cream and a face cloth
4. Paediatric electrodes (use new package)




Storing EEG Cables

All supplies including electrodes, NuPrep, module and leads are in Bay 1 supply room.

After use, slide the two white bars at each end of the electrodes as close to the hub and clips as possible. Let cleaners know to return these bars to this position if they are moved during cleaning.

This will keep the leads from becoming tangled and nurses from losing their minds.

Store on hooks in the South West corner of the Bay 1 supply room (loop the trunk cable on the large hook and the EEG leads on the small hooks).



How to Setting Up Continuous EEG

OBTAIN SUPPLIES

1. EEG module (yellow)
2. ***EEG cable with wires***
3. NuPrep™ cream and a face cloth
4. Paediatric electrodes (new package)



How to Setting Up Continuous EEG

OBTAIN SUPPLIES

1. EEG module (yellow)
2. EEG cable with wires
- 3. *NuPrep™ cream and a face cloth***
4. Paediatric electrodes (new package)



How to Setting Up Continuous EEG

OBTAIN SUPPLIES

1. EEG module (yellow)
2. EEG cable with wires
3. NuPrep™ cream and a face cloth
4. ***Paediatric electrodes (new package)***
repositionable electrodes 2660-3



Prepare Skin


1. Wet a washcloth with warm water (**no soap, cleanser or alcohol**) and scrub skin along hairline and behind ears
2. Rub vigorously with a dry washcloth or towel (to remove residue and roughen epithelial cells)
3. Rub a small amount of Nuprep* into the skin on all areas where electrodes will be applied
4. Remove all traces of Nuprep with a dry washcloth

*Mild abrasive gel that reduces impedance, improves conductivity and reduces artifact



Confirm that the Color Sequence Matches

The black reference electrode is on this side.



Apply Electrodes (diagram next page)

- Apply electrodes just below hair line (sub-hairline); leave enough space to ensure electrode is not in patient's hair which will reduce contact
- Avoid pressing on the centre of the electrode (this can cause a loss of conducting gel)
 1. Place one sub-hairline electrode in the centre of the sub-hairline
 2. Place one electrode in front of each ear
 3. Place one electrode over each temple
 4. Place one electrode between the centre electrode and the temple electrode on each side
 5. Place one electrode behind each ear

Electrode Placement




Identify the CEEG Lead Colors



Connect Electrode to Match Colors as Shown

Right: Even Numbers
Left: Odd numbers



Montage A: Long bipolar (standard lead in critical care)

Monitoring Should Begin Immediately


Monitoring should begin immediately.

If monitoring doesn't begin immediately, make sure that the module is pushed in a connecting to the back of the rack.

If monitoring still doesn't begin, press on the "xylophone" key at the bottom right of the screen an enable the EEG module.



Supported LAN



HR 110/50 **60**

Pulse **60**

PVC 0

ST-I 0.0 ST-aVR 0.0

ST-II -0.1 ST-aVL 0.0

ST-III -0.1 ST-V2 0.2

CVP **(9)**

ICP Mean **(9)**

etCO₂ 50/30 **40**

imCO₂ awRR 30/8 **20**

RR **15**

Delta1 3.5

Theta1 1.5

Alpha1 1.5

Beta1 92.5

Delta2 2.0

Theta2 3.0

Alpha2 2.0

Beta2 93.0

Clock

Touch the change screen and select the EEG waveform. The 4 waves will be displayed as shown here.


20/80 (90)

Man 07:33

SEF1	24.0	SEF3	24.0
MDF1	22.0	MDF3	22.0
PPF1	23.5	PPF3	23.5
SEF2	24.0	SEF4	24.0
MDF2	22.0	MDF4	22.0
PPF2	23.5	PPF4	23.5

7:34

Unsupported LAN



SV Rhythm

HR 110/50 **60**

ART Sys. Mean **120/70 (91)**

Pulse **60**

PVC 0

ST-I 0.0 ST-aVR 0.0
ST-II -0.1 ST-aVL 0.0
ST-III -0.1 ST-V2 0.2

CVP **(9)**

ICP Mean **(9)**

etCO₂ 50/30 **40**

imCO₂ 0
awRR 30/8 **20**

RR **15**

Delta1 2.5
Theta1 3.5
Alpha1 1.5
Beta1 92.5
Delta2 2.0
Theta2 3.0
Alpha2 2.0
Beta2 93.0
Clock

You return to the Dynamic Waves screen if you do not want all waves displayed. You can also choose to change the number of waves continuously displayed. Monitoring continues in the background.

NBP Pulse 60 Man 07:33


Sys. 140
90
Mean **120/80 (90)**
100
60

SEF1	24.0	SEF3	24.0
MDF1	22.0	MDF3	22.0
PPF1	23.5	PPF3	23.5
SEF2	24.0	SEF4	24.0
MDF2	22.0	MDF4	22.0
PPF2	23.5	PPF4	23.5

7:34

From the smart keys at the bottom, select the next page until you find the EEG setup menu.

The monitor automatically detects whether there are 2 (4 montage) or 4 wave (8 montage).



The default montage is “Long Bipolar”.

This is the standard montage for critical care.

The other montages are setup and used by Neuro Critical Care.

EEG Impedance / Montage

Long Bipolar
Long Bipolar
Ipsilat ear
Ref Fpz
Montage D
Montage E

		Impedance
EEG2+	Fpz	2kOhm
EEG2-	A2	2kOhm
EEG3+	F7	2kOhm
EEG3-	T5	2kOhm
EEG4+	F8	2kOhm
EEG4-	T6	2kOhm

TP4

2.54
1.81

15

SEF1	24.0	SEF3 24.0
MDF1	22.0	MDF3 22.0
PPF1	23.5	PPF3 23.5
SEF2	24.0	SEF4 24.0
MDF2	22.0	MDF4 22.0
PPF2	23.5	PPF4 23.5

Delta1	2.5
Theta1	3.5
Alpha1	1.5
Beta1	92.5
Delta2	2.0
Theta2	3.0
Alpha2	2.0
Beta2	93.0
Clock	


7:3

EEG Control

Change Name Change Electrds ImpLimit Down ImpLimit Up

Main Setup Main Screen

to e cycle P



Electrode codes are as follows:

Even numbers right side

Odd numbers left side

F (frontal)

P (parietal)

Fp (frontal polar)

T (temporal)

O (occipital)

A ("zero" or electrode on midline)

C (centre)

EEG1+	Fp1	2kOhm
EEG1-	A1	2kOhm
EEG2+	Fp2	2kOhm
EEG2-	A2	2kOhm
EEG3+	F7	2kOhm
EEG3-	T5	2kOhm
EEG4+	F8	2kOhm
EEG4-	T6	2kOhm

Unsupported LAN



Pulse **60**

PVC 0

ST-I 0.2 ST-aVR -0.1
ST-II 0.0 ST-aVL 0.2
ST-III -0.2 ST-aVF -0.1

CVP **(9)**

ICP Mean **(9)**

etCO₂ **40**

imCO₂ awRR **20**

TP1 **2.34**

TP2 **1.81**

TP3 **2.34**

TP4 **1.81**

SEF1	24.0	SEF3	24.0
MDF1	22.0	MDF3	22.0
PPF1	23.5	PPF3	23.5
SEF2	24.0	SEF4	24.0
MDF2	22.0	MDF4	22.0
PPF2	23.5	PPF4	23.5

Delta1 2.5
Theta1 3.5
Alpha1 1.5
Beta1 92.5
Delta2 2.0
Theta2 3.0
Alpha2 2.0
Beta2 93.0
Clock

EEG Impedance / Montage

Long Bipolar

Long Bipolar

Ipsilat ear


Ref Fpz

Montage D

Montage E

	Impedance
EEG1+	2kOhm
EEG1-	2kOhm
EEG2+	2kOhm
EEG2-	2kOhm
EEG3+	2kOhm
EEG3-	2kOhm
EEG4+	2kOhm
EEG4-	2kOhm

Unsupported LAN



Pulse **60**

PVC 0

ST-I 0.2 ST-aVR -0.1
ST-II 0.0 ST-aVL 0.2
ST-III -0.2 ST-aVF -0.1

CVP **(9)**

ICP Mean **(9)**

etCO₂ **40**

imCO₂ awRR **20**

TP1 **2.34**


TP2 **1.81**

TP3 **2.34**

TP4 **1.81**

SEF1	24.0	SEF3	24.0
MDF1	22.0	MDF3	22.0
PPF1	23.5	PPF3	23.5
SEF2	24.0	SEF4	24.0
MDF2	22.0	MDF4	22.0
PPF2	23.5	PPF4	23.5


Delta1 2.5
Theta1 3.5
Alpha1 1.5
Beta1 92.5
Delta2 2.0
Theta2 3.0
Alpha2 2.0
Beta2 93.0
Clock



EEG Impedance / Montage



Long Bipolar ▾



Impedance Limit:


- EEG1
- EEG2
- EEG3
- EEG4

- 5kOhm
- Fp1-A1
- Fp2-A2
- F7-T5
- F8-T6

Channel	Impedance
Fp1	2kOhm
A1	2kOhm
Fp2	2kOhm
A2	2kOhm
F7	2kOhm
T5	2kOhm
F8	2kOhm
T6	2kOhm

Impedance:

- Measures electrode-to-skin impedance continuously
- Individual electrode impedance is shown in lower table
- Impedance measurements **require that at least 2 electrodes** and a reference electrode must be connected. This electrode relationship is shown by the straight lines in the diagram and in the upper table.



Impedance Code:

- ✓ impedance within range (< 5 kOhm)
- X electrode disconnected
- ↑ impedance above limit
- ? noisy signal

5kOhm
 Fp1-A1
 Fp2-A2
 F7-T5
 F8-T6

Impedance
5kOhm
5kOhm
5kOhm
5kOhm
5kOhm
5kOhm
5kOhm
2kOhm
2kOhm

EEG4+
 EEG4-

F8
 T6

2kOhm
 2kOhm




Checking Impedance:

Use the diagram and table to identify the electrode that requires attention (Example: T6 is > 5).

Re-prep the skin and apply a new electrode. Ensure that the lead is secure within the cable.

If this does not correct the problem, change the “partner electrode”. This is identified by the colored line in the diagram (purple for this example). It can also be identified in the Impedance Limit (the two leads required for EEG channel 4 in this example is F8 and T6).




Long Bipolar

Impedance Limit:	5kOhm
■ EEG1	Fp1-A1
■ EEG2	Fp2-A2
■ EEG3	F7-T5
■ EEG4	F8-T6

Channel	Impedance
EEG1+	Fp1
EEG1-	A1
EEG2+	Fp2
EEG2-	A2
EEG3+	F7
EEG3-	T5
EEG4+	F8
EEG4-	T6


Unsupported LAN



You can return to the Dynamic Waves screen if you do not want EEG continuously displayed. You can also choose to change the number of waves being displayed. Monitoring will still continue in the background and the waves will be available in full disclosure from the central station.



Unsupported LAN



To stop displaying any wave, touch on the waveform and choose "Change Wave".

EEG Wave

- Freeze Wave
- Wave Scale : 100µV
- Setup EEG
- Change Speed
- Change Wave**

SV Rhythm HR 110/50 **60**

Pulse **60**

PVC 0

ST-I 0.2 ST-aVR -0.1
ST-II 0.0 ST-aVL 0.2
ST-III -0.2 ST-V6 0.4

CVP **(9)**

ICP Mean **(9)**

TP1 **2.34**

TP2 **1.81**

TP3 **2.34**

TP4 **1.81**

etCO2 **40**

imCO2 awRR **20**

RR **15**

Delta1 2.5
Theta1 3.5
Alpha1 1.5
Beta1 92.5
Delta2 2.0
Theta2 3.0
Alpha2 2.0
Beta2 93.0
Clock

NBP Pulse 60 Man 07:33 NBP mmHg 07:33 **120/80(90)**

Sys. 140
90
Mean 100
60

SEF1	24.0	SEF3	24.0
MDF1	22.0	MDF3	22.0
PPF1	23.5	PPF3	23.5
SEF2	24.0	SEF4	24.0
MDF2	22.0	MDF4	22.0
PPF2	23.5	PPF4	23.5

7:36



Unsupported LAN



DEMO

SV Rhythm

HR

110
50

60

Pulse



60

PVC 0

ST-I 0.3 ST-aVR -0.2
ST-II 0.0 ST-aVL 0.3
ST-III -0.3 ST-V4 1.4

CVP



(9)

ICP

Mean

10

0

(9)

etCO₂

50

30

40

imCO₂

awRR

30

8

20

RR




15

Delta1 2.5
Theta1 3.5
Alpha1 1.5
Beta1 92.5
Delta2 2.0
Theta2 3.0
Alpha2 2.0
Beta2 93.0
Clock

To stop displaying any wave, touch on the waveform and choose "Change Wave".

Change Waves

- Wave 1 : Any RT Wave
- Wave 2 : Any RT Wave
- Wave 3 : Any EEG
- Wave 4 : Any EEG
- Wave 5 : Any EEG
- Wave 6 : Any EEG



TP1 2.34
TP2 1.81
TP3 2.34
TP4 1.81



NBP Pulse 60
Sys. 140
90
Mean 100
60
120/80 (90)

Man 07:33 NBP 07:33 120/80(90)

SEF1	24.0	SEF3	24.0
MDF1	22.0	MDF3	22.0
PPF1	23.5	PPF3	23.5
SEF2	24.0	SEF4	24.0
MDF2	22.0	MDF4	22.0
PPF2	23.5	PPF4	23.5

7:37

✓ Unsupported LAN



Choose the wave you want to remove. "Wave 6" is the 6th waveform on the screen starting with the ECG. In this example, the 6th wave is EEG3.

NBP Pulse 60
 Sys. 140
 90
 Mean 100
 60
120/80 (90)


07:33 120/80(90) SEF1 24.0 SEF3 24.0
 MDF1 22.0 MDF3 22.0
 PPF1 23.5 PPF3 23.5
 SEF2 24.0 SEF4 24.0
 MDF2 22.0 MDF4 22.0
 PPF2 23.5 PPF4 23.5

TP2 1.81
 TP3 2.34
 TP4 1.81

Delta1 2.5
 Theta1 3.5
 Alpha1 1.5
 Beta1 92.5
 Delta2 2.0
 Theta2 3.0
 Alpha2 2.0
 Beta2 93.0
 Clock

7:39

✓ Unsupported LAN



Arrow up or down until you see "Blank". This will turn the display off. Monitoring of this wave will continue in the background.

Change Waves

Wave 1	: Any RT Wave	Resp
Wave 2	: Any RT Wave	EEG1
Wave 3	: Any EEG	EEG2
Wave 4	: Any EEG	EEG3
Wave 5	: Any EEG	EEG4
Wave 6	: Any EEG	Blank
		Primary Lead
		SecondaryLead

NBP Pulse 60
 Sys. 140
 90
 Mean 100
 60
120/80 (90)

07:33 120/80(90) SEF1 24.0 SEF3 24.0
 MDF1 22.0 MDF3 22.0
 PPF1 23.5 PPF3 23.5
 SEF2 24.0 SEF4 24.0
 MDF2 22.0 MDF4 22.0
 PPF2 23.5 PPF4 23.5

Delta1 2.5
 Theta1 3.5
 Alpha1 1.5
 Beta1 92.5
 Delta2 2.0
 Theta2 3.0
 Alpha2 2.0
 Beta2 93.0
 Clock

7:39

✓ Unsupported LAN




If you want to return to the full 4 wave display, go to Change Screen and choose EEG.

Notice the EEG* This is a modified version of the EEG screen. For example, if you had turned 2 of the 4 EEG waves off, this would be your modified screen.

You can go back and forth between the default EEG screen (EEG) and the modified screen (EEG*).



Unsupported LAN



HR	110/50	60	Pulse	60
ART Sys.	120/70	(91)	PVC	0
ART Mean			ST-I	0.0
			ST-II	-0.1
			ST-III	-0.1
			ST-aVR	0.0
			ST-aVL	0.0
			ST-V2	0.2
			CVP	(9)
PAP Sys.	28/15	(21)	ICP Mean	(9)
PAP 40			etCO ₂	40
PAP 20			imCO ₂	0
TP1	2.34		awRR	20
TP2	1.81		RR	15
TP3	2.34		Delta1	2.5
TP4	1.81		Theta1	3.5
			Alpha1	1.5
			Beta1	92.5
			Delta2	2.0
			Theta2	3.0
			Alpha2	2.0
			Beta2	93.0
			Clock	7:34

NBP Pulse 60 Man 07:33

Sys. 140
90
Mean 100
60

120/80 (90)

SEF1	24.0	SEF3	24.0
MDF1	22.0	MDF3	22.0
PPF1	23.5	PPF3	23.5
SEF2	24.0	SEF4	24.0
MDF2	22.0	MDF4	22.0
PPF2	23.5	PPF4	23.5