

# Alpha-Iduronidase, Leukocytes

---

## Additional Names

MPSI, Hurler Syndrome, Scheie Syndrome

## Orderable

AIDUL

## Specimen Requirements

**Specimen Types:** Whole Blood

**Containers:**

- 2 x 6mL Dark Green top (Sodium Heparin, no separator gel) Vacutainer

**Minimum Testing Volume:** 8 – 10mL

## Specimen Handling

**Collection Instructions:** Sample must be collected so that it will be received in lab no later than 3:00 pm of collection day.

## Storage and Transportation

Transport specimens to the laboratory immediately after collection. Keep cool, but not frozen.

## Referred in Specimens Only

A normal control sample can be added by the sender

## Testing Location

Biochemical Genetics Lab

## Testing Schedule

As required

## Method of Analysis

Plate Method



# Reference Intervals

Test	Range
Leukocyte	5.5 – 15.3 nmol/hr/mg protein

