

AMIODARONE /
DESETHYLAMIODARONE
SERUM

Orderable - AMIO

Turnaround Time: 10 days

STAT: N/A

Specimen:

6 mL Red top Vacutainer tube

Collection Information:

Minimum volume required is 0.5 mL.

Reference Ranges:

Amiodarone	Desethylamiodarone
2.30 -3.90 µmol/L	2.30-3.90 µmol/L

Storage and Shipment:

Separate and freeze. Ship frozen and do not allow time to thaw.



Laboratory:
Toxicology/Special
Chemistry Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
HPLC-MS/MS



Test Schedule:
Monday – Friday
0800-1600