

AMPHETAMINE
SCREEN,
URINE

Orderable - AMPSU

Turn Around Time: 24 hours

STAT: 1 hour

Alternate Name(s):

Amphetamines
Methamphetamine
Dextroamphetamine
MDMA
Ecstasy



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Kinetic Interaction of
Microparticles in Solution
(KIMS)



Test Schedule:
As required

Specimen:

Random urine

Collection Information:

Collect random urine in a sterile container.

Reference Ranges:

Negative

Interpretive Comments:

A positive result for the Amphetamine Screen indicates that any of the following drugs may be present in urine alone or in combination: methamphetamine, amphetamine, MDMA (Ecstasy) or dextroamphetamine, at a minimum total concentration of 1000 ng/ml. This result is to be regarded as preliminary/presumptive unless subsequently verified by a definitive confirmation test (available upon request by calling the Toxicology Laboratory).