

ANTI GA1,
IGM,
SERUM

Orderable - AGA1

Turn Around Time: 50 days

Alternate Name(s):

AGA1
Anti-Asialo IgM
Anti Asialo GM1 (GA1)
Anti-Asialo GM1
Anti-GA1 IgM
ASGM1
IgM vs GA1
IgM vs Asialo-GM1 (GA1)
IgM antibodies to Asialo GM1

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube or 6 mL Red top Vacutainer tube	0-2 years: 0.5 mL Red top Microtainer 2-10 years: 0.5 mL Red top Microtainer

Collection Information:

N.B: This test is an out-of-country request. Pre-approval by the Immunology section head (Dr. L. Yang ext. 35768), or Biochemist on-call in her absence Monday-Friday 0800-1600 hrs is required prior to collection.

Minimum serum volume required is 1 mL in each of 2 aliquots.

Reference Ranges:

Titre: < 3000 Units

Interpretive Comments:

Associated with chronic motor neuropathy syndromes.



Laboratory:
Core Lab



Requisition:
[TESTING LAB](#)
[REQUISITION FORM](#)



Method of Analysis:
ELISA



Test Schedule:
Referred out monthly



Referred Out Location:
Washington University
Medical School
660 S. Euclid, Box 8111
St. Louis, Missouri
63110
USA



ANTI GA1,
IGM,
SERUM

Comments:

This test is available exclusively to SJH/LHSC physicians.

For more information, please contact:

Senior Technologist
Immunology Lab
(519) 685-8500 ext. 35541

Storage and Shipment:

Refrigerate up to 2 days, otherwise freeze at -20°C for longer storage.