

ANTI TR, SERUM

Orderable - IMMIX

Turn Around Time: 50 days

Alternate Name(s):

Purkinje Cell Cytoplasmic Antibody Type Tr

Specimen:

5 mL Gold or 6 mL Red top Vacutainer tube

Collection Information:

Minimum serum volume is 1.0 mL in each of 2 aliquots.

Note: This test is available exclusively to LHSC/SJHC physicians.
For more information, please contact the Senior Technologist in the Immunology Lab.
(519)-685-8500 ext. 35541

Reference Ranges:

Negative

Interpretive Comments:

The presence of ant-Tr is usually associated with paraneoplastic cerebellar degeneration and Hodgkin's disease.

Comments:

This test is available exclusively to LHSC/SJHC physicians.
For more information, please contact:
Senior Technologist
Immunology Lab, LHSC
(519)-685-8500 ext. 35541



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
IFA



Test Schedule:
Referred out monthly



Referred Out Location:
Washington University
Medical School
660 S. Euclid, Box 8111
St. Louis, Missouri
63110
USA



ANTI
TR,
SERUM

Storage and Shipment:

Refrigerate up to 2 days, otherwise freeze at
-20°C for longer storage.