

## ANTIENDOMYSIAL ANTIBODIES

# Orderable - EMA

Turnaround Time: 25 days

### Alternate Name(s):

EMA  
Anti Mycelial

### Specimen:

6 mL Red top Vacutainer tube

### Collection Information:

A minimum volume of 1 mL of serum is required for testing.

### **Note:**

This test is not to be used for celiac disease screening unless the patient is TTGAb negative and IgA deficient. **This test is approved for Dr. Andreanne Zizzo, Dr. Tejas Desai and Dr. Dhandapani Ashok.** All others will require biochemist approval to order this test; approval can be requested by calling the clinical biochemist at extension 77626 or the biochemist on-call.

### Reference Ranges:

Negative

### Specimen Stability:

Refrigerated: 7 days  
Frozen: >7 days



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Immunofluorescence  
Assay (IFA)



**Test Schedule:**  
Referred out Monday-  
Thursday



**Referred Out Location:**  
[In-Common Laboratories](#)