



BENZODIAZEPINES,
URINE

Orderable - BENSU

Turn Around Time: 24 hours

STAT: 1 hour

Alternate Name(s):

Benzo Screen

Specimen:

Random urine

Collection Information:

Collect random urine in a sterile container.

Reference Ranges:

Negative

Interpretive Comments:

A positive result for the Benzodiazepine Screen indicates that the drug may be present in urine at a minimum concentration of 200 ng/mL, and is to be regarded as preliminary/presumptive unless subsequently verified by a definitive confirmation test (available upon request by calling the Toxicology Laboratory).



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Kinetic interaction of
Particles in Solution
(KIMS)



Test Schedule:
As required