

BILIRUBIN - DIRECT

Orderable - BILD

Turn Around Time: 4 hours

STAT: 1 hour

Alternate Name(s):

Direct
Conjugated

Specimen:

Adult	Pediatric
4.5 mL Green top Vacutainer	0-2 years: 0.5 mL Green top Microtainer 2-10 years: 2 mL Green top Vacutainer tube

Collection Information:

For pediatric samples, a minimum of 0.2 mL is required.

Turn off phototherapy lights prior to drawing sample.

Reference Ranges:

Age	Reference Interval ($\mu\text{mol/L}$)
0-<15 days	≤ 9
15 days-<1 year	≤ 4
1-<9 years	≤ 2
9-<13 years	≤ 4
13-<19 years	≤ 5
Adult	≤ 3



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Timed endpoint Diazo
method



Test Schedule:
As required