





Pathology and Laboratory Medicine

# Orderable - BIU / BIU24

Turnaround Time: 10 days

#### **Specimen:**

24-hour urine collected in an unused 24-hour urine container or random urine

## **Collection Information:**

**Trace Elements Collection Requirements** 

### **Reference Ranges:**

SI Units:

Random Urine: 0-32 nmol/mol creatinine

24 Hour Urine: 0-0.48 nmol/d

#### **Conventional Units:**

Random Urine: 0-59 ng/g creatinine

24 Hour Urine 0-0.100 μg/d

Reference Ranges are based on Non-Occupationally exposed population.

# **Interpretive Comments:**

Find Interpretive Comment and Clinical Information here:

#### **Specimen Stability:**

Ambient: 14 days

Refrigerated: 11 months Frozen: 11 months



**Laboratory:** Trace Elements Lab

**Requisition:** TRACE ELEMENTS REQUISITION



Method of Analysis: ICP-MS/MS



**Test Schedule:** Batched analysis