

Pathology and Laboratory Medicine

CADMIUM, URINE

Orderable - CADU / CADU24

Turnaround Time: 10 days

Specimen:

24-hour urine collected in an unused 24-hour urine container or random urine

Collection Information:

0

Laboratory: Trace Elements Lab

Trace Elements Collection Requirements



Requisition: TRACE ELEMENTS REQUISITION

Method of Analysis:

ICP-MS/MS

Reference Ranges:

<u>SI Units:</u>

Random Urine: 0-1.4 μmol/mol creatinine

24 Hour Urine: 0-17.8 nmol/d

•_• ⊞

Test Schedule: Batched analysis

Conventional Units:

Random Urine: 0-1.4 μg/g creatinine

24 Hour Urine: 0-2.0 μg/d

Interpretive Comments:

Find Interpretive Comment and Clinical Information here:

Specimen Stability:

Ambient: 14 days Refrigerated: 11 months Frozen: 11 months