

CAFFEINE,
SERUM

Orderable - CAFF

Turnaround Time: 5 days

Specimen:

6 mL Red top Vacutainer tube
or
4.5 mL Green (Lithium Heparin) top Vacutainer tube

Collection Information:

Minimum volume of serum or plasma required is 0.5 mL.

Spin and aliquot serum or plasma and freeze immediately and send to the lab.

Reference Ranges:

25-80 $\mu\text{mol/L}$

Critical Information Required:

Please specify if the sample is from an infant.



Laboratory:
Toxicology/Special
Chemistry Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
HPLC-MS/MS



Test Schedule:
Monday – Friday
0800-1600