

VARICELLA
ZOSTER
SEROLOGY
IgG/IgM

Orderable – VZVG/VZVM

Turnaround Time: Up to 5 days

Specimen Type and Collection Information:

Please review the Public Health Ontario Laboratory Test Information Index:

Varicella - Serology

Specimen Type	Minimum Volume	Collection Kit
Whole Blood or Serum	5.0 mL blood or 1.0 mL serum	Blood, clotted – Gold vacutainer tube (SST)

Storage and Shipment:

The specimen must be refrigerated.
The specimen must be received within 24 hours of collection.



Laboratory:
Microbiology



Requisition:
See electronic order



Method of Analysis:
VZVG – Multiplex Flow
Immunoassay
VZVM - CLIA



Test Schedule:
Referred out weekdays



Referred to:
Public Health Ontario
Laboratory
1200 Commissioners Rd E
London, Ontario
519-455-9310