

CLOBAZAM
/ DESMETHYL
CLOBAZAM,
SERUM

Orderable - CLB

Turnaround Time: 5 days

Specimen:

6 mL Red top Vacutainer tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Minimum volume of serum required for testing is 0.5 mL.

Reference Ranges:

Clobazam	Desmethyl Clobazam
No therapeutic range has been established, however levels of: 0.15-1.00 µmol/L has been suggested	2.8-14.0 µmol/L has been suggested

Storage and Shipment:

Separate and freeze. Ship frozen



Laboratory:
Toxicology/Special
Chemistry Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
HPLC-MS/MS



Test Schedule:
Monday – Friday
0800-1600