

CREATININE
URINE -
RANDOM

Orderable - CREU

Turn Around Time: 4 hours

Specimen:

Random urine

Collection Information:

Collect random urine in a sterile container with no preservative.

Reference Ranges:

Random Urine (1st morning):

Male	Female
3.5-25 mmol/L	2.6-20 mmol/L

Interpretive Comments:

Creatinine excretion is included as part of a creatinine clearance test or may be ordered separately.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Enzymatic method



Test Schedule:
Monday - Friday
0800 - 1600 h