



DIHYDROPTERIDINE
REDUCTASE,
DRIED
BLOOD
SPOT

Orderable - DHPR

Alternate Name(s):

PKU-DHPR Deficient

Specimen:

Filter Paper with 6 mL Green (Sodium Heparinized) top Vacutainer

Collection Information:

Saturate a 1.5 X 4 cm strip of a PKU-type card with whole blood (Na-Heparin) *or* heel prick blood; air dry thoroughly.

Reference Ranges:

7-22 U/g Hb

Storage and Shipment:

Sample must be sent to lab **immediately**.



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)



Test Schedule:
As required