



FRUCTOSAMINE,  
SERUM

## Orderable - FRUCTO

Turnaround Time: 5 days

---

### Specimen:

6 mL Red top Vacutainer tube

### Collection Information:

Minimum volume required is 1 mL

### Reference Ranges:

205 - 285  $\mu\text{mol/L}$

### Storage and Shipment:

Store and send sample frozen.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Nitroblue



**Test Schedule:**  
Referred out Monday-  
Thursday



**Referred Out Location:**  
[In-Common Laboratories](#)