

GALACTOSE-1-
PHOSPHATE
URILTRANSFERASE -
QUALITATIVE,
DRIED BLOOD
SPOT

Orderable - G1PUDBS

Alternate Name(s):

Galactosemia Screen

Specimen:

Dried blood spots on a Newborn Screening-type card.

Blood can be from a heel prick, finger prick, or Green top (Sodium or Lithium Heparinized) Vacutainer tube.

Collection Information:

On a Newborn Screening-type card, saturate each circle completely (i.e blood dot should be clearly visible on both front and back of card).

Air dry horizontally at room temperature. Do **NOT** heat the card to dry!

Once the blood spots are completely dry (approximately 3 hours) seal the card in a plastic bag and send to Biochemical Genetics at room temperature.

Reference Ranges:

See report



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)



Test Schedule:
As required