

GLUTETHIMIDE
SERUM

Orderable – GENX: then specify test

Turn Around Time: 5 days

Alternate Name(s):

Doriden

Specimen:

5 mL Gold top Vacutainer tube

Collection Information:

Avoid gel-separator tubes

Minimum volume required is 1 mL.

Reference Ranges:

Therapeutic
0.9 - 3.7 µmol/L
Toxic
Greater than 45 µmol/L

Storage and Shipment:

Store and send cold.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
LC-MS/MS



Test Schedule:
Referred out Monday -
Thursday



Referred Out Location:
[In-Common Laboratories](#)