



HEPATITIS B SURFACE ANTIGEN

Orderable - HpBsAg

Turn Around Time: Up to 48 hours

Alternate Name(s):

HBsAg
HpBsAg

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 0.5 mL Gold top Microtainer 2-10 years: 2 mL Gold top tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

See report

Interpretive Comments:

Earliest serum indicator of acute Hepatitis B infection. Persistent elevation may indicate Hepatitis B carrier state.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Electrochemiluminescence
Immunoassay



Test Schedule:
24/7