

 HLA
B27

Orderable - HLA B27

Turnaround Time: 14 days

Specimen:

2 x 4mL Lavender (EDTA) top Vacutainer tubes

Collection Information:

Samples that are being shipped into the laboratory can be drawn Monday to Friday.

Reference Ranges:

See report

Interpretive Comments:

HLA B27 Positive indicates that the patient expresses HLA B27 antigen on his/her cells.

HLA B27 Negative indicates that the patient does not express HLA B27 antigen on his/her cells.

Storage and Shipment:

Keep blood at room temperature and deliver within 72 hours of collection.



Laboratory:
Transplant Lab



Requisition:
[TRANSPLANT
LABORATORY
REQUISITION](#)



Method of Analysis:
PCR-SSP



Test Schedule:
As required Monday to
Friday except Stat
holidays