

HLA CELIAC INVESTIGATION

Orderable - HLA CELIAC

Turnaround Time: 14 days

Specimen:

2 x 4 mL Lavender (EDTA) top Vacutainer tubes

Collection Information:

Samples that are being shipped into the laboratory can be drawn Monday to Friday.

Reference Ranges:

See report

Interpretive Comments:

HLA DQ2 (DQA1*05) or DQ8 Positive/Negative

Celiac is associated with DQ8 or with DQ2 in combination with DQA1*05

Storage and Shipment:

Keep blood at ambient temperature and deliver within 72 hours.



Laboratory:
Transplant Lab



Requisition:
[TRANSPLANT
LABORATORY
REQUISITION](#)



Method of Analysis:
HLA Class I and Class II
 typings are performed
 using molecular
 techniques (RSSO or SSP)



Test Schedule:
As required Monday to
Friday except Stat
holidays