

JAPANESE B ENCEPHALITIS VIRUS SEROLOGY

Orderable - ARBOSER

Turnaround Time: 28 days

Specimen and Collection Information:

Please review the Public Health Ontario Laboratory Test Information Index:

[Japanese B Encephalitis Virus](#)

Specimen Type	Minimum Volume	Collection Kit
Serum	1.0 mL serum	Gold Vacutainer tube (SST)

Critical Information Required:

Ensure the Arbovirus Information Intake form includes the following:

- Onset date
- Relevant travel history
- History of Japanese encephalitis virus vaccination or yellow fever vaccination

Specify **Japanese B Encephalitis** in the **Specimen Description** field when placing order in PowerChart and on the Arbovirus Information Intake form .

Storage and Shipment:

The specimen must be refrigerated.
The specimen must be received within 24 hours of collection.



Laboratory:
Microbiology



Requisition:
[Arbovirus \(Non-Zika\)
Information Intake Form](#)
This requisition **MUST**
accompany the specimen



Method of Analysis:
ELISA
PRNT



Test Schedule:
Referred out weekdays



Referred to:
Public Health Ontario
Laboratory
1200 Commissioners Rd E
London, Ontario
519-455-9310