



Pathology and Laboratory Medicine

LEAD, URINE

Orderable - PBU / PBU24

Turnaround Time: 10 days

Specimen:

24-hour urine collected in new 24-hour urine container or random urine

Collection Information:

Trace Elements Collection Requirements

Reference Ranges:

SI Units:

Random Urine:

0-1.0 μmol/mol creatinine

24 Hour Urine: 0-19.3 nmol/d

Conventional Units:

Random Urine: $0-1.8 \mu g/g$ creatinine

24 Hour Urine: 0-4.0 μg/d

Reference Ranges are based on Non-Occupationally exposed population.

Interpretive Comments:

Find Interpretive Comment and Clinical Information here:

Specimen Stability:

Ambient: 14 days

Refrigerated: 11 months Frozen: 11 months



Laboratory: Trace Elements Lab

<u>≡</u>

Requisition: TRACE ELEMENTS REQUISITION





Test Schedule:Batched analysis