



LEFLUNOMIDE

Orderable - LEFLU

Turn Around Time: 3 weeks

Alternate Name(s):

Arava
Teriflunomide

Specimen:

6 mL Red top or 4 mL Lavender (EDTA) top Vacutainer tube

Collection Information:

Serum separator gel tubes are not acceptable.

Freeze plasma/serum aliquots prior to shipment.

Reference Ranges:

50-100 ug/mL

Note: For couples planning a pregnancy, confirmation of teriflunomide concentrations (either father or mother receiving therapy) of <0.02 ug/mL in two separate tests at least 14 days apart is highly recommended.

Comments:

If the results are to be faxed back from the referral laboratory, please indicate so clearly on the request.

Storage and Shipment:

Separate and freeze. Ship on ice packs or frozen.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
HPLC



Test Schedule:
Weekly, test is performed
on Thursdays



Referred Out Location:
CHUM – Laboratoire
Regroupe
1000 rue St-Denis
Montreal, Quebec
H2X 0C1
Telephone: 514-890-8318
Fax: 514-412-7420