

LIDOCAINE SERUM/ PLASMA

Orderable - LIDO

Turnaround Time: 7 days

Alternate Name:

Xylocaine

Specimen:

Adult	Pediatric
6 mL Red top Vacutainer tube	0-2 years: 0.5 mL Red top Microtainer 2-10 years: 3 mL Red top Vacutainer tube
Plasma (heparin or EDTA) is also acceptable	

Collection Information:

Blood collection tubes with separator gels are not acceptable for this test.

Collect trough specimen within 1 hour prior to next dose.

Minimum volume of serum (or plasma) required is 200 µL.

Reference Ranges:

1.5-5.0 mg/L

Storage and Shipment:

Ambient: 5 days
Refrigerated: 7 days
Frozen: 30 days

Store and ship serum or plasma frozen.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Immunoassay



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)