

LIPASE

Orderable - LIP

Turn Around Time: 24 hours

Specimen:

Adult	Pediatric
4.5 mL Green (Lithium Heparin) top Vacutainer	0-2 years: 0.5 mL Green Microtainer 2-10 years: 2 mL Green top tube
Serum from a 5 mL Gold top or 6 mL Red top Vacutainer tube is also acceptable	



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Enzymatic colorimetric



Test Schedule:
As required

Reference Ranges:

13-60 U/L

Interpretive Comments:

Elevated in acute pancreatitis. Not very specific for pancreatic disorders. Also increased in primary biliary cirrhosis, chronic renal failure and with hemodialysis.