

LIPID PROFILE

Orderable - LIPIDS

Turnaround Time: 24 hours

Alternate Name(s):

Cholesterol, Triglyceride, HDL-C, LDL-C

Specimen:

| Adult | Pediatric |
|---|--|
| 4.5 mL Green (Lithium Heparin) top Vacutainer tube | 0-2 years: 0.5 mL Green Microtainer 2-10 years: 1 mL Green top tube |
| Red, Gold, or Lavender (EDTA) top tubes are also acceptable | |

Collection Information:

12-14 hour fast recommended. Non-fasting specimens acceptable.

Reference Ranges:

See individual tests for reference ranges.

Interpretive Comments:

See individual tests for interpretive comments.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
See individual tests:
Cholesterol, Triglycerides,
HDL-C & LDL-C



Test Schedule:
As required