

MAGNESIUM,  
WHOLE BLOOD

## Orderable - MGBLD

Turnaround Time: 10 days

### Specimen:

6 mL K2-EDTA Royal Blue Vacutainer tube

### Collection Information:

[Trace Elements Collection Requirements](#)

### Reference Ranges:

SI Units (Reported on Patient Chart):	Conventional Units:
1.04-1.49 mmol/L	25.3-36.3 mg/L

Reference Ranges are based on Non-Occupationally exposed population.

### Interpretive Comments:

Find Interpretive Comment and Clinical Information [here](#):

### Specimen Stability:

Ambient: 22 days  
Refrigerated: 14 months  
Frozen: 14 months



**Laboratory:**  
Trace Elements Lab



**Requisition:**  
TRACE ELEMENTS  
REQUISITION



**Method of Analysis:**  
ICP-MS/MS



**Test Schedule:**  
Batched analysis