

MUMPS
SEROLOGY
IgG/IgM

Orderable – MUMG/MUMM

Turnaround Time: Up to 5 days

Please review the Public Health Ontario Laboratory Services Test Information Index:

[Mumps Serology](#)

Specimen Type and Collection Information:

Specimen Type	Minimum Volume	Collection Kit
Blood or Serum	5.0 mL blood or 1.0 mL serum	Blood, clotted – Gold vacutainer tube (SST)

Acute: Collect the acute sample within 7 days after the onset of the rash.

Convalescent: Collect the convalescent sample 7 to 10 days after the acute

Storage and Shipment:

Specimen must be refrigerated.

Specimen must be received within 24 hours of collection.



Laboratory:
Microbiology



Requisition:
See electronic order



Method of Analysis:
IgG – Multiplex Flow
Immunoassay
IgM - CLIA



Test Schedule:
Referred out weekdays



Referred to:
Public Health Laboratory
1200 Commissioners Rd E
London, Ontario
519-455-9310