

# NITRAZEPAM SERUM

## Orderable - NITRAZ

Turnaround Time: 5 days

### Alternate Name(s):

Mogadon

### Specimen:

2 x 6 mL Red top Vacutainer tube

or

2 x 4.5 mL Green (Lithium Heparin) top Vacutainer tube

### Collection Information:

Collect blood aseptically in a Vacutainer tube.

Minimum volume of serum or plasma required is 3 mL.

Spin and aliquot serum or plasma and freeze immediately and send to the lab.

### Reference Ranges:

100-800 nmol/L

### Storage and Shipment:

Separate and freeze. Ship frozen and do not allow time to thaw



**Laboratory:**  
Toxicology/Special  
Chemistry Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
LC-MS/MS



**Test Schedule:**  
Monday--Friday  
0800-1600