

PHOSPHOLIPASE
A2 RECEPTOR
(PLA2R)
ANTIBODIES

Orderable – PLA2R

Turnaround Time: 14 days

Alternate Name(s):

Anti-PLA2R

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube or 6 mL Red top Vacutainer tube	0-2 years: 0.5 mL Red or Gold top Microtainer 2-10 years: 2 mL Red top Vacutainer tube

Collection Information:

Testing restricted to nephrologists.
Minimum volume of serum is 0.5 mL.

Reference Ranges:

Negative:	<14 RU/mL
Borderline:	14 - 19 RU/mL
Positive:	>19 RU/mL

Interpretive Comments:

Anti-PLA2R antibodies can be detected in the serum of up to 70-80% of patients with primary membranous nephropathy.

Comments:

This test is available exclusively to SJHC/LHSC nephrologists.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION




Method of Analysis:
ELISA



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)



PHOSPHOLIPASE
A2 RECEPTOR
(PLA2R)
ANTIBODIES

Storage and Shipment:

Store and send frozen.

Specimen Stability:

Refrigerated: 14 days

Ambient: Unacceptable