

POTASSIUM
URINE -
RANDOM

Orderable - KU

Turn Around Time: 4 hours

STAT: 1 hour

Specimen:

Random urine

Collection Information:

Collect a minimum of 5 mL of random urine in a sterile container.

Reference Ranges:

30-120 mmol/L

Comments:

RBC hemolysis may cause interference with results.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Indirect potentiometry



Test Schedule:
As required