

SALICYLATE,
SERUM

Orderable - SAL

Turn Around Time: 4 hours

STAT: 1 hour

Alternate Name(s):

Salicylic Acid
ASA

Specimen:

Adult	Pediatric
4.5 mL Green top Vacutainer tube	0-2 years: 0.5 mL Green top Microtainer 2-10 years: 2 mL Green top tube
Serum from a 5 mL Gold top or 6 mL Red top Vacutainer tube is also acceptable	

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

Analgesic/Antipyretic:
0.2-0.7 mmol/L
Anti-inflammatory:
1.1-2.2 mmol/L
Toxic:
>2.3 mmol/L
Critical Value:
>2.3 mmol/L



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Colorimetric



Test Schedule:
As required