

# SODIUM, 24-HOUR URINE

## Orderable - NAU24

Turn Around Time: 4 hours

STAT: 1 hour

### Specimen:

24-hour urine collection

### Collection Information:

#### [Collection of 24-Hour Urine Sample](#)

No preservative is required.

Keep timed specimen container in refrigerator or on ice during the timed period.

Please complete the laboratory label attached to the collection bottle to include patient's name, collection start time and date and collection end time and date.

### Reference Ranges:

40-217 mmol/d

### Storage and Shipment:

Store timed specimen container in refrigerator or on ice during the timed period.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Indirect potentiometry  
using glass sodium  
electrodes.



**Test Schedule:**  
Weekdays