

SULFATIDES  
URINE

## Orderable - SULFTU

Turn Around Time: 1 month

### Alternate Name(s):

Metachromatic Leukodystrophy  
MLD

### Specimen:

Random urine

### Collection Information:

Collect a minimum of 4 mL of random urine in sterile, screw top container and send frozen to lab.

Keep at -20°C.

### Reference Ranges:

0-170 µg/L

### Storage and Shipment:

Ship sample frozen on dry ice.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Quantification of  
sulfatides of varying fatty  
acid chain length by LC-  
MS/MS.



**Test Schedule:**  
Referred out as required



**Referred Out Location:**  
[The Hospital for Sick  
Children](#)  
Rapid Response  
Laboratory  
170 Elizabeth Street,  
Room 3642  
Toronto, ON  
M5G 2G3  
Canada  
Phone: 1-855-381-3212