





Pathology and Laboratory Medicine

# Orderable - SULFU / SULFU24

Turnaround Time: 10 days

#### **Specimen:**

24-hour urine collected in new 24-hour urine container or random urine

## **Collection Information:**

**Trace Elements Collection Requirements** 

### **Reference Ranges:**

SI Units:

Random Urine: 1.1-3.6 mmol/mol creatinine

24 Hour Urine: 7.8-25.0 mmol/d

#### **Conventional Units:**

Random Urine: 324-1011 mg/g creatinine

24 Hour Urine: 250-802 mg/d

Reference Ranges are based on Non-Occupationally exposed population.

# **Interpretive Comments:**

Find Interpretive Comment and Clinical Information <a href="here">here</a>:

#### **Specimen Stability:**

Ambient: 14 days

Refrigerated: 11 months Frozen: 11 months



**Laboratory:** Trace Elements Lab

**≡** ×-

Requisition: TRACE ELEMENTS REQUISITION



Method of Analysis: ICP-MS/MS



**Test Schedule:**Batched analysis