

# SYSTEMIC SCLEROSIS PROFILE

## Orderable - .ANANP

Turnaround Time: 5 days

### Alternate Name(s):

Anti Scl-70  
Anti CENP A  
Anti CENP B  
Anti RP11  
Anti RP155  
Anti Fibrillarin  
Anti NOR90  
Anti Th/To  
Anti PM-Scl100  
Anti PM-Scl75  
Anti Ku  
Anti PDGFR  
Anti Ro-52

### Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube or 6 mL Red top Vacutainer tube	0-2 years: Red 0.5 Microtainer 2-10 years: 2 mL Red top
<b>EDTA, Heparin or Citrate plasma are also acceptable</b>	

### Collection Information:

Minimum volume required for testing is 0.5 mL

### Reference Ranges:

Negative

### Interpretive Comments:

Immunoblot results will be reported as either negative or positive for specific antibodies as the following:

- (+) Borderline
- 1+ Weak Positive
- 2+ Positive
- 3+ Strong Positive



**Laboratory:**  
Clinical Immunology Lab




**Requisition:**  
[CLINICAL IMMUNOLOGY  
REQUISITION](#)



**Method of Analysis:**  
Immunoblot



**Test Schedule:**  
Batch analysis



# SYSTEMIC SCLEROSIS PROFILE

## Comments:

**This test is available exclusively to Rheumatologists and Respiriologists at LHSC/SJHC and accepted from referred in locations.**

## Specimen Stability:

Sample can be stored for up to 14 days at 4°C. For longer storage, store at -20°C.