

TRYPTASE,  
SERUM/  
PLASMA

## Orderable - TRY

Turn Around Time: 15 days

### Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube or 6 mL Red top Vacutainer tube	0-2 years: Red 0.5 mL Microtainer 2-10 years: 2 mL Red top tube
Plasma is also acceptable. 4.5 mL Green top Vacutainer tube	



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Fluorescent Enzyme  
Immunoassay (FEIA)



**Test Schedule:**  
Referred out Monday-  
Thursday



**Referred Out Location:**  
[In-Common Laboratories](#)

### Collection Information:

Minimum volume required for testing is 1 mL.

Separate ASAP.

To assess anaphylaxis, collect specimen between 15-180 minutes after suspected anaphylactic event. To assess systemic mastocytosis or mast cell activation syndrome, collect the specimen at any time.

### Reference Ranges:

3.8 – 11.4 µg/L

### Storage and Shipment:

Store and ship frozen. Avoid repeated freeze/thaw cycles.