

VANADIUM, URINE

Orderable – VU / VU24

Turnaround Time: 10 days

Specimen:

24-hour urine collected in new 24-hour urine container or random urine

Collection Information:

Trace Elements Collection Requirements



Laboratory:
Trace Elements Lab



Requisition:
TRACE ELEMENTS
REQUISITION



Method of Analysis:
ICP-MS/MS



Test Schedule:
Batched analysis

Reference Ranges:

SI Units:

Random Urine:
0-176 nmol/mol creatinine

24 Hour Urine:
0.39-3.93 nmol/d

Conventional Units:

Random Urine:
0-79 µg/g creatinine

24 Hour Urine:
0.020-0.200 µg/d

Reference Ranges are based on Non-Occupationally exposed population.

Interpretive Comments:

Find Interpretive Comment and Clinical Information [here](#):

Specimen Stability:

Ambient: 14 days
Refrigerated: 11 months
Frozen: 11 months