



WHOLE  
BLOOD TOXIC  
PANEL,  
PLASMA

## Orderable - TETOXBLD (Toxic Panel)

Turnaround Time: 10 days

---

### Specimen:

6 mL K2 EDTA Royal Blue Vacutainer tube

### Collection Information:

[Trace Elements Collection Requirements](#)

### Reference Ranges:

[Find individual Reference Ranges here:](#)

Reference Ranges are based on Non-Occupationally exposed population.

### Interpretive Comments:

Find Interpretive Comment and Clinical Information [here:](#)



**Laboratory:**  
Trace Elements Lab



**Requisition:**  
TRACE ELEMENTS  
REQUISITION



**Method of Analysis:**  
ICP-MS/MS



**Test Schedule:**  
Batched analysis