ZINC PROTOPORPHYRIN, WHOLE BLOOD

Orderable - ZPP

Turnaround Time: 10 days

Specimen:

6 mL K2 EDTA Royal Blue Vacutainer tube

Collection Information:

Trace Elements Collection Requirements



Requisition: TRACE ELEMENTS REQUISITION

Laboratory: Trace Elements Lab

Reference Ranges:

μmol/mol heme 0-70

Reference Ranges are based on Non-Occupationally exposed population.

Interpretive Comments:

Find Interpretive Comment and Clinical Information here:



Test Schedule: Batched analysis

Method of Analysis: Proto Fluor London Health Sciences Centre

Pathology and Laboratory Medicine