

MAPROTILINE,
SERUM/
PLASMA

Orderable - MAP

Turnaround Time: 5 days

Alternate Name(s):

Ludiomil

Specimen:

6 mL Red top Vacutainer tube

Lavender or Pink (EDTA) top tubes are also acceptable

Avoid blood collection tubes with separator gels (e.g., Gold top).

Light Green (Li-Heparin) top tubes are also NOT acceptable.

Collection Information:

Minimum volume of serum or EDTA plasma required is 3 mL.

To monitor therapy, draw trough specimen prior to next dose or 10 - 12 h after last drug administration. Separate as soon as possible.

Reference Ranges:

Please note change in units, therapeutic range, and toxic range as of October 18, 2022.

Maprotiline	75-130 ng/mL
-------------	--------------

Critical Value:

Toxic:

Maprotiline	> 220 ng/mL
-------------	-------------

Storage and Shipment:

Separate serum as soon as possible.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
LC-MS/MS



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)