



NEUTRALIZING
ANTIBODIES
TO
INTERFERON-
BETA

Orderable - NAB

Turn Around Time: 20 days

Alternate Name(s):

NAB

Specimen:

5 mL Gold top Vacutainer tube

Collection Information:

Test only available to Dr. M. Kremenchutsky, Dr. C. Hyson and Dr. Sarah Morrow of the M.S. Clinic at University Hospital, London Health Sciences Centre.

Special Processing:

Centrifuge at 4 °C, transfer serum to cryovials and store at 4 °C.

Critical Information Required:

The minimum information required to process this test is the patients OHIP #, Name and D.O.B.

Storage and Shipment:

Store serum at 4 °C and ship with ice packs.



Laboratory:
Core Lab



Requisition:
[UBC Neuro-Immunology Laboratory Requisition](#) must accompany the specimen and is to be filled in by the MS clinic



Method of Analysis:
Reporter Gene Assay



Test Schedule:
Referred out Mondays



Referred Out Location:
Faculty of Medicine
Koerner Pavilion
Room S-157 (2nd floor),
2211 Wesbrook Mall
Vancouver, BC Canada
V6T 2B5
Tel (Office) 604 822 7696
| (Lab) 604 822 7175
Fax 604 822 0758
Website:
<https://bcneuro.ca/>
Email:
neuroimm.lab@ubc.ca