

POTASSIUM  
24-HOUR  
URINE

## Orderable – KU24

Turn Around Time: 4 hours

STAT: 1 hour

### Specimen:

24-hour urine

### Collection Information:

#### Collection of 24-Hour Urine Sample

No preservative required in 24-hour urine container.

Please complete the laboratory label attached to the collection bottle to include patient's name, collection start time and date and collection end time and date.

### Reference Ranges:

30-130 mmol/d

### Comments:

RBC hemolysis may cause interference with results.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Indirect potentiometry



**Test Schedule:**  
Weekdays