

UREA, 24- HOUR URINE

Orderable – UU24

Turn Around Time: 4 hours

Specimen:

24-hour urine

Collection Information:

[Collection of 24-Hour Urine Sample](#)

No preservative necessary, but specimen should be kept refrigerated during collection.

Please complete the laboratory label attached to the collection bottle to include patient's name, collection start time and date and collection end time and date

Reference Ranges:

170-580 mmol/24 hours

Storage and Shipment:

Keep specimen refrigerated then send to the lab on ice packs.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Enzymatic Photometric



Test Schedule:
Weekdays