

ABO TITRE

Orderable - A, B TITRE IGM, IGG

Turn Around Time: 4 hours

STAT: 2 hours

Alternate Name(s):

Alloantibody
Antibody Titre

Specimen:

6 mL Pink (EDTA) top Vacutainer tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

See report

Interpretive Comments:

ABO antibody titres are done to measure the amount of anti-A and/or anti-B in major ABO incompatible bone marrow transplant recipients. Anti-A titres are done for group O kidney transplant recipients if a group A2 kidney donor is going to be considered.

Critical Information Required:

Diagnosis

Storage and Shipment:

Transport to the Blood Transfusion Laboratory within 24 hours of collection at room temperature.



Laboratory:
Blood Transfusion Lab



Requisition:
[BLOOD TRANSFUSION REQUISITION](#) or
Electronic order



Method of Analysis:
Anti-Human Globulin test



Test Schedule:
Daily
Urgent, if indicated