



DIRECT ANTIGLOBULIN TEST

Orderable - DIRECT ANTIGLOBULIN TEST

Turn Around Time: 45 minutes

STAT: 15 minutes

Alternate Name(s):

DAT
Coombs' Test

Specimen:

6 mL Pink (EDTA) top Vacutainer tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

See report

Interpretive Comments:

Usually ordered only when special hematological work-up is being done.

Storage and Shipment:

Transport to the Blood Transfusion Laboratory within 24 hours of collection at room temperature.



Laboratory:
Blood Transfusion
Laboratory



Requisition:
[BLOOD TRANSFUSION
REQUISITION](#) or
Electronic order



Method of Analysis:
Direct antiglobulin test



Test Schedule:
Daily