

TISSUE FOR
IMMUNOFLUORESCENCE
(IF) – SKIN,
CONJUNCTIVAL
BIOPSY, BUCCAL
MUCOSA

Orderable - Surgical Pathology Request

Turnaround Time: <2 weeks

STAT: 6 hours

Specimen:

Tissue biopsy

Collection Information:

Contact Pathology Lab at ext. 36307 prior to collection for Immunofluorescent fixative.

Place biopsy immediately in 5 mL vial with Immunofluorescent fixative.

Reference Ranges:

Negative

Interpretive Comments:

Lupus Band Test
Pemphigus and Bullous pemphigoid antibodies



Laboratory:
Pathology Lab



Requisition:
PowerChart: E-order
choosing appropriate
specimen.
See [Identification of
Clinical Specimens](#)



Method of Analysis:
Direct
immunofluorescence



Test Schedule:
Tuesday and Thursday