

MYASTHENIA
GRAVIS
PANEL
(anti-AChR
and MuSK)

Orderable - ACHR/MUSK

Turnaround Time: 7 days

Alternate Name(s):

ACHR Antibodies
AChRab
MuSK Antibodies
MuSK Ab



Laboratory:
Clinical Immunology Lab



Requisition:
[CLINICAL IMMUNOLOGY
REQUISITION](#)



Method of Analysis:
Cell-based assay (CBA)



Test Schedule:
Batch analysis

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube or 6 mL Red top Vacutainer tube	0-2 years: 0.5 mL Red top Microtainer 2-10 years: 2 mL Red top Vacutainer tube
EDTA, heparin or citrated plasma are also acceptable	

Collection Information:

Minimum volume of serum is 1 mL (2 mL preferred)

Reference Ranges:

Negative

Interpretive Comments:

Antibodies associated with Myasthenia Gravis.

Comments:

This test is restricted to neurologists and ophthalmologists at LHSC/SJHC.
Accepted from referred in locations.

MYASTHENIA
GRAVIS
PANEL
(anti-AChR
and MuSK)

Sample Stability:

Refrigerated: 14 days

Frozen: 30 days